



DENDRIDIAG® BF

**OPERATING MODE FOR QUANTIFICATION OF TOTAL FLORA ON
SURFACES BY ATP-METRIE**

- LUMINOMÈTER KIKKOMAN PD30 -

VERSION: V2018-10

CONTENTS

GL BIOCONTROL overview	3
What is ATP-metry?	4
Why use ATP-metry for microbiological monitoring?	5
Equipment needed	6
Protocol for quantification of total flora on a surface	8
Result analysis	10
Calculation of biomass quantity	10
Data management with Excel	10
Result interpretation	12
Troubleshooting	13
Controls	14
F.A.Q.	15
Contact	21

GL BIOCONTROL overview

GL BIOCONTROL specializes in **environmental risk management** and has an expertise in sanitary engineering along with biological monitoring of water and surfaces.

Our clients are environment professionals: industrialists, industry operators, water treatment companies, laboratories, study design engineer and public authorities.

We have several skills including: studies, research and development, analysis, product development and professional training. Through these areas of expertise, GL BIOCONTROL:

- **Develops** risk management tools (ATP-metry kits for total flora quantification, DNA extraction purification kits, real time PCR amplification kits, electropositive membranes for viruses...).
- **Uses** methodologies and innovative tools to study the microbial world (qPCR, NGS, ATP-metry...).
- **Studies** ecosystems to anticipate and prevent public health risks, in particular linked to *Legionella* and *Pseudomonas* genus (risk assessment, microbiology diagnostics, ATP cartography...).
- **Advises** water sector professionals on how to manage their facilities in order to reduce public health risks as well as improve the environmental footprint (water, treatment products and energy conservation).
- **Trains** environmental professionals on microbiological risk management and laboratory techniques.

GL BIOCONTROL offers everything you need to quantify total flora in waters, on surfaces and in the air by **ATP-metry**: measurement kits **DENDRIDIAG®**.




What is ATP-metry?

Adenosine triphosphate (ATP) is the major intermediary energy required in most cellular metabolism reactions. Every living cell produces and consumes ATP. This coenzyme, specific to living environments, proves the existence of living organisms.

Quantifying ATP equates to **quantifying total microorganisms** (or total flora).

To perform this type of assay, the light emitted by the enzymatic reaction of **bioluminescence** using luciferin and firefly luciferase is measured (see below).

4



ATP, in the presence of a luciferin/luciferase complex with a catalyst, releases energy in the form of light. By measuring the amount of light emitted using a **luminometer**, we deduce the **quantity of ATP in picogram per square centimeter**. The total flora, expressed in equivalent bacteria per square centimeter, is calculated from the following:

$$1 \text{ picogram} \approx 1\,000 \text{ bacteria.}$$

The ATP-metry measurement method is a **field test** whose result is obtained in few minutes.

Why use ATP-metry for microbiological monitoring?

Contamination or degradation of surface microbial quality is caused by one or all of the following: dubious cleaning and disinfection procedures, nature of the biocide used, low quality of the surface material. There is no definitive solution to eradicate microbiological issues. Only an approach starting with the **observation** of the surfaces state, followed by the **monitoring** of the corrective action effects can result in effective management of microbiological risks.

Most of the regulatory texts relating to microbiological risk management of surfaces, require that the facility manager uses **indicators** to follow and anticipate a microbiological shift of their facilities, in order to avoid contaminants such as (*Campylobacter*, *Pseudomonas* or *Salmonella* for example).

The monitoring indicator should be a technology that is **rapid**, **reliable**, **easy to use** and **economic** in order to frequently analyze the facility.


Among the different microbiological indicators, the most frequent are heterotrophic plate count at 22°C or 36°C, quantitative PCR, qualitative ATP-metry and **quantitative ATP-metry**.

Quantitative ATP-metry is one of the best indicators for biological monitoring. Using quantitative ATP-metry, you will:

- **Control the facility cleanliness:** manage biofouling, improve health risk management, Avoid production shutdown due to non-compliance with the regulatory requirement (*Pseudomonas*, *Campylobacter*...).
- **Assess operating procedure efficiency:** validate efficiency of cleaning (bio-dispersant) and disinfection (biocide) procedures, validate treatment strategy.
- **Manage biofilm formation in a pipe:** control and assess the fouling state of the pipe network, start and validate cleaning and disinfection procedures.
- **Identify the critical areas:** Determine critical areas with strong biomass growth, adapt cleaning and disinfection strategy, highlight malfunctions.

Equipment needed

Kit of reagents for 100 measurements



Product	Quantity	
Dropper bottle DENDRIDIAG® BF	5	
Dropper bottle EXTRACTANT	5	
Dropper bottle STANDARD 1000	2	






Reagents (**DENDRIDIAG® BF**, **EXTRACTANT** and **STANDARD 1000**) should be stored in the dark in a freezer (-18°C). In this way, they can be kept for at least 12 months. After first use, the reagents will be preferentially refrozen. Or else, they can be kept refrigerated (between 3 and 8°C) for a maximum of 8 consecutive weeks.

6

Kit of consumables for 100 measurements

Product	Quantity	
Sterile swabs	100	
Disposable sterile polypropylene test tubes	100	

Equipment

Product	Quantity	
Stainless steel sampling template (20 cm ²)	1	
Luminometer KIKKOMAN PD30 or equivalent	1	
Tube holder for PD30	1	
Fridge (3 to 8°C) *	1	
Freezer (≈ -18°C) *	1	

* For good conservation of the reagents, it is necessary to store them in a freezer or at least in a fridge.

Protocol for quantification of total flora on surfaces

Phase 1: installation

1. Thaw a dropper bottle of each reagent: **DENDRIDIAG® BF**, **EXTRACTANT** and **STANDARD 1000**. Bring them to room temperature (above 18°C),
2. Prepare the plastic consumables (swab and test tube) and the sampling template,
3. Turn on the luminometer and wait 10 seconds for the device calibration.






8



*To perform the analysis, the reagents **DENDRIDIAG® BF**, **EXTRACTANT** and **STANDARD 1000** must be at room temperature (between 18°C and 25°C) to ensure a maximal enzyme efficiency.*

Phase 2: sampling

1. Make sure the reagent **EXTRACTANT** is close to room temperature (> 18°C) and put 8 drops of **EXTRACTANT** in the middle of the test tube,
2. Get the swab out of its packaging being careful not to touch the cotton tip with your fingers or with the lab bench,
3. Plunge the cotton tip in the test tube containing the **EXTRACTANT**, turn the swab in the test tube to humidify it,
4. Place the sampling template on the surface you want to control,
5. Scrub the surface defined by the sampling template (20cm²) with the swab. Pass at least two times on the surface on two different directions,






From this step, no break time is allowed.

Phase 3: measurement



Make sure the reagent is close to room temperature (> 18°C). Warm up the reagent in your hand if necessary.

6. Add 2 drops of **DENDRIDIAG® BF** in the middle of the test tube containing the **EXTRACTANT**,
7. Put the swab in the test tube and homogenize the mix by turning the swab,
8. Break the swab while maintaining it in the test tube. Do not touch the cotton tip of the swab with your fingers,
9. Fix the tube to the tube holder,
10. Place them in the luminometer and press the ENTER button to start measurement,
11. After 10 seconds of measurement, write down the R1 result in RLU (Relative Light Unit),
12. Immediately, get the tube holder with the test tube out of the luminometer,
13. Take the test tube off the tube holder and add one drop of **STANDARD 1000** in the middle of the test tube. The dropper bottle must not touch the tube while adding the drop,
14. Correctly homogenize the mix by turning the swab inside the test tube to get a good standardization,
15. Fix the tube to the tube holder,
16. Place them in the luminometer and press the ENTER button to start measurement,
17. After 10 seconds of measurement, write down the R2 result in RLU (Relative Light Unit).



Result analysis

Calculation of total flora quantity

The intracellular ATP concentration is expressed in picogram ATP per square centimeter. To obtain the result, execute the following operations:

Calculation of the standard (in RLU/pgATP):

$$\text{STANDARD} = \frac{R2 - R1}{1\ 000}$$

Calculation of the biomass value (in pgATP/cm²):

$$[\text{ATP}] = \frac{R1}{\text{STANDARD} \times S}$$

10

With:

R1 = result obtained on the sample in RLU,

R2 = result obtained on the sample + **STANDARD 1000** in RLU,

STANDARD = value of the **STANDARD 1000** in RLU/pgATP,

S = surface sampled (defined at 20 cm² with the sampling template).

It is possible to convert the ATP concentration (in pgATP/cm²) to equivalent bacteria per square centimeter (eq. bact./cm²) using the following rule: 1 picogramme ATP ≈ 1 000 bacteria.

For example, 5 pg/cm² ≈ 5000 eq.bact./cm².

Data management with Excel

An Excel file is supplied by GL BIOCONTROL to automatically perform the calculations presented above. The file must be used to monitor your facility and interpret the values obtained.

After each analysis, sampling date or location, surface sampled, R1 and R2 values measured by the luminometer must be entered in the table supplied. Only the grey columns must be completed.

The table automatically performs the calculation. Results are given in:

- **picogram** per square centimeter (pg ATP/cm²),
- total flora in **equivalent bacteria** per square centimeter (eq.bact./ cm²),
- total flora in **logarithm** per square centimeter (LOG eq.bact./cm²).


By default, we use the result in logarithm equivalent bacteria per square centimeter. However, you can use one of the three units.

A		B		C		D		E		F		G		H		I		J		K	
1		2		3		4		5		6		7		8		9		10		11	
12		13		14		15		16		17		18		19		20		21		22	
23		24		25		26		27		28		29		30		31		32		33	
34		35		36		37		38		39		40		41		42		43		44	
45		46		47		48		49		50		51		52		53		54		55	
56		57		58		59		60		61		62		63		64		65		66	
67		68		69		70		71		72		73		74		75		76		77	
78		79		80		81		82		83		84		85		86		87		88	
89		90		91		92		93		94		95		96		97		98		99	
100		101		102		103		104		105		106		107		108		109		110	
111		112		113		114		115		116		117		118		119		120		121	
122		123		124		125		126		127		128		129		130		131		132	
133		134		135		136		137		138		139		140		141		142		143	
144		145		146		147		148		149		150		151		152		153		154	
155		156		157		158		159		160		161		162		163		164		165	
166		167		168		169		170		171		172		173		174		175		176	
177		178		179		180		181		182		183		184		185		186		187	
188		189		190		191		192		193		194		195		196		197		198	
199		200		201		202		203		204		205		206		207		208		209	
210		211		212		213		214		215		216		217		218		219		220	
221		222		223		224		225		226		227		228		229		230		231	
232		233		234		235		236		237		238		239		240		241		242	
243		244		245		246		247		248		249		250		251		252		253	
254		255		256		257		258		259		260		261		262		263		264	
265		266		267		268		269		270		271		272		273		274		275	
276		277		278		279		280		281		282		283		284		285		286	
287		288		289		290		291		292		293		294		295		296		297	
298		299		300		301		302		303		304		305		306		307		308	
309		310		311		312		313		314		315		316		317		318		319	
320		321		322		323		324		325		326		327		328		329		330	
331		332		333		334		335		336		337		338		339		340		341	
342		343		344		345		346		347		348		349		350		351		352	
353		354		355		356		357		358		359		360		361		362		363	
364		365		366		367		368		369		370		371		372		373		374	
375		376		377		378		379		380		381		382		383		384		385	
386		387		388		389		390		391		392		393		394		395		396	
397		398		399		400		401		402		403		404		405		406		407	
408		409		410		411		412		413		414		415		416		417		418	
419		420		421		422		423		424		425		426		427		428		429	
430		431		432		433		434		435		436		437		438		439		440	
441		442		443		444		445		446		447		448		449		450		451	
452		453		454		455		456		457		458		459		460		461		462	
463		464		465		466		467		468		469		470		471		472		473	
474		475		476		477		478		479		480		481		482		483		484	
485		486		487		488		489		490		491		492		493		494		495	
496		497		498		499		500		501		502		503		504		505		506	
507		508		509		510		511		512		513		514		515		516		517	
518		519		520		521		522		523		524		525		526		527		528	
529		530		531		532		533		534		535		536		537		538		539	
540		541		542		543		544		545		546		547		548		549		550	
551		552		553		554		555		556		557		558		559		560		561	
562		563		564		565		566		567		568		569		570		571		572	
573		574		575		576		577		578		579		580		581		582		583	
584		585		586		587		588		589		590		591		592		593		594	
595		596		597		598		599		600		601		602		603		604		605	
606		607		608		609		610		611		612		613		614		615		616	
617		618		619		620		621		622		623		624		625		626		627	
628		629		630		631		632		633		634		635		636		637		638	
639		640		641		642		643		644		645		646		647		648		649	
650		651		652		653		654		655		656		657		658		659		660	
661		662		663		664		665		666		667		668		669		670		671	
672		673		674		675		676		677		678		679		680		681		682	
683		684		685		686		687		688		689		690		691		692		693	
694		695		696		697		698		699		700		701		702		703		704	
705		706		707		708		709		710		711		712		713		714		715	
716		717		718		719		720		721		722		723		724		725		726	
727		728		729		730		731		732		733		734		735		736		737	
738		739		740		741		742		743		744		745		746		747		748	
749		750		751		752		753		754		755		756		757		758		759	
760		761		762		763		764		765		766		767		768		769		770	
771		772		773		774		775		776		777		778		779		780		781	
782		783		784		785		786		787		788		789		790		791		792	
793		794		795		796		797		798		799		800		801		802		803	
804		805		806		807		808		809		810		811		812		813		814	
815		816		817		818		819		820		821		822		823		824		825	
826		827		828		829		830		831		832		833		834		835		836	
837		838		839		840		841		842		843		844		845		846		847	
848		849		850		851		852		853		854		855		856		857		858	
859		860		861		862		863		864		865		866		867		868		869	
870		871		872		873		874		875		876		877		878		879		880	
881		882		883		884		885		886		887		888		889		890		891	
892		893		894		895		896		897		898		899		900		901		902	
903		904		905		906		907		908		909		910		911		912		913	
914		915		916		917		918		919		920		921		922		923		924	
925		926		927		928		929		930		931		932		933		934		935	
936		937		938		939		940		941		942		943		944		945		946	
947		948		949		950		951		952		953		954		955		956		957	
958		959		960		961		962		963		964		965		966		967		968	
969		9																			


Result interpretation

Warning and alarm thresholds were established based on our experience. These thresholds should be refined based on the first results obtained on your facilities. The following arrows will help you **interpret the results**:


Operating swimming pools (in LOG/cm²):



Surfaces of food processing industry (in pgATP/cm²):



Surfaces of water pipes (in pgATP/cm²):



We consider that when a measurement is:

- **Below the warning threshold**, the surface is under control,
- **Between the warning and the alarm threshold**, the surface does not present an immediate biohazard. A corrective action is recommended if 3 consecutive measurements are above the warning threshold,
- **Above the alarm threshold**, the surface is not under control. A quick corrective action is recommended.

In the Excel file, the result is displayed in green, orange or red depending if it is under the warning threshold, between the warning and the alarm threshold or above the alarm threshold.

In order to help understanding the results, we encourage you to write down any relevant information in the "user's comment" column.

Troubleshooting

The Excel file automatically detects 3 types of error and displays a message explaining the issue identified. When an error is highlighted, please refer to the troubleshooting table below.

If a measurement must be repeated after an error, we advise you to replace the values on the same line to avoid errors on the graph sheet.

Problem	Possible cause and correction
“Low sensitivity. If necessary, sample a larger surface.” displayed in the Excel file.	<p>The reagent DENDRIDIAG® BF is not sufficiently active (out-of-date, degraded or too cold) to obtain high sensitivity.</p> <p>Warm up the reagents DENDRIDIAG® BF and EXTRANCTANT and sample a larger surface if necessary. If the problem remains, perform a Control of the reagent efficiency (cf. page 14).</p>
“Control the mixing of the STANDARD 1000, the temperature and condition of the reagents.” displayed in the Excel file.	<p>Standardization was not successful. Tap the bottom of the test tube on a flat surface and homogenize the mix by turning the swab in the test tube. Restart the measurement.</p> <p>If the problem remains, the reagent DENDRIDIAG® BF is not sufficiently active (out-of-date, degraded or too cold). Warm up the reagents DENDRIDIAG® BF and EXTRANCTANT and perform a Control of the reagent efficiency (cf. page 14).</p>
“Sample highly contaminated. If necessary, sample a smaller surface.” displayed in the Excel file.	<p>The ATP concentration of the sample is too high. Restart the analysis by sampling a smaller surface.</p>

If you have troubles with our kit, do not hesitate to contact us by phone or email:

Yannick FOURNIER
Sales engineer
9, avenue de l'Europe - Cap Alpha
34 830 CLAPIERS (FRANCE)
Phone: +33 (0)6 33 64 42 29
Phone: +33 (0)9 67 39 35 20
Email: y.fournier@gl-biocontrol.com

Clément FAYE
Research engineer
9, avenue de l'Europe - Cap Alpha
34 830 CLAPIERS (FRANCE)
Phone: +33 (0)6 72 70 46 98
Phone: +33 (0)9 67 39 35 20
Email: c.faye@gl-biocontrol.com

Controls

Control of the luminometer contamination

a) Test:

- Fix an empty test tube to the tube holder,
- Place it in the luminometer and press the ENTER button,
- The result should be less or equal to 2 RLU.

b) Protocol to be followed in case of contamination:

With a cotton swab, wipe the internal surfaces of the measurement chamber.

14

Control of the reagent contamination

a) Test:

- In a test tube, put 2 drops of DENDRIDIAG® BF and 4 drops of EXTRACTANT,
- Fix the test tube to the tube holder,
- Place it in the luminometer and press the ENTER button,
- The result should be less or equal to 5 RLU.

b) Protocol to be followed in case of contamination:

Discard the contaminated reagents and select a new bottle of DENDRIDIAG® BF and of EXTRACTANT.

Control of the reagent efficiency

a) Test:

- In a test tube, put 2 drops of DENDRIDIAG® BF, 4 drops of EXTRACTANT and 1 drop of STANDARD 1000 (reagent temperature must be above 18°C),
- Fix the test tube to the tube holder,
- Properly homogenize the tube,
- Place it in the luminometer and press the ENTER button,
- For a good efficacy of the reagents, the result should be higher than 130 RLU.

b) Protocol to be followed in case of degradation:

Discard the reagents and select a new bottle of DENDRIDIAG® BF and of EXTRACTANT.

F.A.Q.

GENERAL POINTS ON ATP AND DENDRIDIAG® KITS

What is ATP?

Adenosine triphosphate (ATP) is the major intermediary energy required in most cellular metabolism reactions. Every living cell produces and consumes ATP. This coenzyme, specific to living environments, proves the existence of living organisms.

What do I measure with ATP-metry?

In water, quantifying ATP equates to quantifying total microorganisms (or total flora). ATP-metry is a biomolecular technique, based on bioluminescence. The measurement is done using a luminometer.

What do I measure with the DENDRIDIAG® kits?

With DENDRIDIAG® kits and filtration, only intracellular ATP is measured. It corresponds to the ATP found inside the living cells representative of living bacteria.

Extracellular ATP is also found in sample as a free molecule in the sample. It comes from dead or dying microorganisms. Filtration eliminates free ATP. Without filtration, total ATP is measured: intracellular and extracellular ATP.

Which microorganisms are lysed by the DENDRIDIAG® reagent?

DENDRIDIAG® kits preferentially lyses bacteria, cyanobacteria and amoeba. For total lysis of all microorganisms (fungi, yeast, algae...) consult GL BIOCONTROL.

At what temperature should I use the reagents?

To ensure a maximal enzyme efficiency, DENDRIDIAG® reagent and STANDARD 1000 must be used at room temperature (18°C - 25°C).

How and how long can I store the reagents?

All ATP-metry dropper bottles (DENDRIDIAG® reagents, STANDARD 1000 and EXTRACTANT) must be stored in the dark in a freezer (-18°C). In this way, they can be kept for at least 12 months. After first use, the reagents will be preferentially refrozen or, if not, kept refrigerated (between 3 and 8°C) for 8 consecutive weeks.

Stored at room temperature, the reagents are stable less than one week.

How and how long can I keep the plastic consumables?

Plastic consumables must be stored in a dry area at room temperature. Their expiration date is displayed on their individual packaging.

I forgot the DENDRIDIAG® reagent at room temperature. What should I do?

For good stability of the kit, all the reagents must be stored in a freezer (-18°C). If you forgot the reagents at room temperature for a few days, you can perform a control of the reagent efficiency. To do so, refer to the paragraph *Control of the reagent efficiency* page 14.

16

In which areas of application can use the DENDRIDIAG® kits?

Industrial Water (IW): cooling system, process water circuit, water supply system for industrial purposes...

Sanitary Water (SW): drinking water supply system, water network of spa facilities...

Ultra-Pure Water (UPW): water loop system for medical, pharmaceutical or microelectronics use, water networks under microbiological control...

Surface (BF): food processing, cooling tower, water supply system, pools...

Air (AIR): aeraulic network, hospitals, offices, high risks industries like composting, methanation, farming...

OPERATING MODE

What surface dimension should I sample?

Our kit is supplied with a sampling template which defines a surface of 20 cm². With this sampling template, you get a good reproducibility and a significant sampling of the surface studied which leads to reliable results. Nevertheless, if you want to increase the measurement sensitivity, you can perform two sampling of the same surface by moving the sampling template. Always write down the surface sampled.

17

LUMINOMETER

My tube holder for PD30 luminometer is dirty. How to clean it?

An excess of foam or a bubble in the upper part of the tube can dirty the tube holder. In this case, you can use a paper towel impregnated with alcohol 70% of water to clean the tube holder.

I forgot to write down the R1 and R2 values.

It is possible to retrieve the results in RLU in the luminometer. To do so, turn on the luminometer and after the 10 seconds of calibration, press the up arrow to get the last values obtained.

Luminometer PD30 rings and does not complete calibration.

PD30 calibration is done when the luminometer is empty and the cap correctly closed. If the luminometer rings, make sure the tube holder is out and the cap is well closed.

Luminometer C110 displays « OVERSCALE ».

Luminometer C110 is very sensitive. The message "OVERSCALE" is displayed when the sample is strongly contaminated and the luminometer cannot measure the RLU. If you analyze liquids, restart the measurement by filtrating 1/10th of the volume.

RESULT INTERPRETATION

Which unit should I use to express my results?

Picogram ATP per square centimeter (pgATP/cm²) is the real unit.

For a better understanding of the results, it is possible to use the unit equivalent bacteria per square centimeter (eq.bact./cm²) using the scientific consensus 1 pgATP ≈ 1 000 bacteria. This result is not rigorously true because the ATP concentration varies from microorganism to microorganism and also differs following the metabolic state of the bacteria.

18

Why express my results in LOG?

Results are often expressed in LOG to simplify the interpretation. Indeed, we consider that two results are significantly different if there is a difference of 1 LOG between the two values.

The following table shows you the conversion of eq.bact./cm² in LOG:


eq.bact./cm ²	LOG(eq.bact./cm ²)
10	1
100	2
1 000	3
10 000	4
100 000	5

The Excel table automatically gives you this value.


What reference limits should I use for my surface analysis?

Based on our experiment, we established the following warning and alarm thresholds:


Operating swimming pools (in LOG/cm²):



Surfaces of food processing industry (in pgATP/cm²):



Surfaces of water pipes (in pgATP/cm²):




19

These thresholds should be refined based on the first results obtained on your network.

Why is the R2 value inferior to the R1 value?

Signal loss was estimated to 6% to 8% per minute. If measurement takes a long time and there are several minutes between ATP extraction and reading of R2, standardization will happen during the signal degrowth phase. It is often observed on sample highly contaminated.

We advise you to restart the analysis by sampling a smaller surface.



Why the message “Low sensitivity. If necessary, sample a larger surface.” Is displayed in the Excel file?

The reagent DENDRIDIAG® BF is not sufficiently active (out-of-date, degraded or too cold) to obtain high sensitivity.

Warm up the reagents DENDRIDIAG® BF and EXTRANCTANT and sample a larger surface if necessary. If the problem remains, perform a *Control of the reagent efficiency* (cf. page 14).

Why the message “Control the mixing of the STANDARD 1000, the temperature and condition of the reagents.” Is displayed in the Excel file?

Standardization was not successful. Tap the bottom of the test tube on a flat surface and homogenize the mix by turning the swab in the test tube. Restart the measurement.

If the problem remains, the reagent DENDRIDIAG® BF is not sufficiently active (out-of-date, degraded or too cold). Warm up the reagents DENDRIDIAG® BF and EXTRANCTANT, and perform a *Control of the reagent efficiency* (cf. page 14).

Why the message “Sample highly contaminated. If necessary, sample a smaller surface.” is displayed in the Excel file.

The ATP concentration of the sample is too high. Restart the analysis by sampling a smaller surface.

Contact

For further information of assistance on interpretation of results, on the protocol or for commercial information, contact by email or by phone:

Yannick FOURNIER

Sales engineer

9, avenue de l'Europe - Cap Alpha

34 830 CLAPIERS (FRANCE)

Phone: +33 (0)6 33 64 42 29

Phone: +33 (0)9 67 39 35 20

Email: y.fournier@gl-biocontrol.com

Clément FAYE

Research engineer

9, avenue de l'Europe - Cap Alpha

34 830 CLAPIERS (FRANCE)

Phone: +33 (0)6 72 70 46 98

Phone: +33 (0)9 67 39 35 20

Email: c.faye@gl-biocontrol.com

21

A tutorial video of the protocol is available on the USB key supplied with the luminometer or on our website in the tab Products – ATP-metry kit for surfaces:

www.gl-biocontrol.com

4 easy ways to order

@ by email at contact@gl-biocontrol.com,

by fax at + 33 (0)9 55 25 40 31,

by phone at + 33 (0)9 67 39 35 20,

by mail at GL BIOCONTROL - 9, avenue de l'Europe, Cap Alpha - 34 830 CLAPIERS (FRANCE).

>>> www.gl-biocontrol.com

