

CONTROLE THE QUALITY OF YOUR ANALYZES

ISO 10705-2 requires the analysis of a positive control for each series of sample.

The Φ X174 Positive Control is a freeze-dried water-soluble biological material, containing a specified number of Φ X174 phages.

THE **ΦX174** POSITIVE CONTROL FOR A MAXIMUM EFFICIENCY



Save time thanks to ready-to-use tablets



Comply with the standard by using a positive control in each series of analysis



Only one step simplified by an instant solubilization



Benefit from high precision material

PRACTICAL INFORMATION



without defrosting
 rehydrate with sterile water



Storage
- temperature: -20°C
- shelf life: 9 months



Content 4 bags of 5 units



Target Around 100 PFU per vial

NEED INFORMATION? CONTACT US!



+ 33 (0)9 67 39 35 20



contact@gl-biocontrol.com



RAPID MICROBIAL DIAGNOSTIC



GL BIOCONTROL 5 avenue de l'Europe - Hélioparc -34 830 CLAPIERS (FRANCE)

www.gl-biocontrol.com