



DENDRIDIAG® UPW

ATP tests for ultrapure water



MEASURE THE QUANTITY OF BACTERIA
IN YOUR WATER IN LESS THAN 10 MINUTES

TAKE CONTROL

of your water circuit.
Avoid production defects and
the resulting costs of poor
quality.

ANTICIPATE

microbiological deviations
in water quality with a kit
offering fast results and high
sensitivity.

SAVE

on manpower and treatment
products. Trigger disinfection
operations in a reasoned
manner.

MASTER YOUR WATER NETWORK



Monitor

the microbial quality of your
ultra-pure water loops.



Conduct

your interventions
intelligently.



ASSESS

the efficacy of preventive and
corrective treatments.



Identify

the critical points
in no time at all.

THEY TRUST US



WHAT IS DENDRIDIAG®?



A turnkey solution to quickly determine the quantity of bacteria in your ultra-pure water networks.


▲ Watch the video

- Dialysis water
- Water for medical and pharmaceutical preparation.
- Water for microelectronic production.




Get the result in a second thanks to the **web app!**

A FAST AND RELIABLE ANALYSIS



INDIVIDUAL SUPPORT




KIT CONTENT



The reagents

- 3 DENDRIDIAG droppers
- 3 EXTRACTANT droppers
- 1 STANDARD dropper



The consumables

- 60 extension tubes and pipettes
- 60 50 ml syringes and 0.45 µm filters
- 60 test tubes



The measuring device

- 1 luminometer KIKKOMAN C110
- 1 microtube rack
- 1 USB flash drive

NEED INFORMATION? CONTACT US



Head office

5, avenue de l'Europe - Hélioparc
34 830 CLAPIERS (FRANCE)

Grand Est Agency

8, rue Alfred Kastler
Maison de l'entreprise
54 320 MAXÉVILLE (FRANCE)



+ 33 (0)9 67 39 35 20



contact@gl-biocontrol.com



Join us on **Linkedin**

GLBIOCONTROL

RAPID MICROBIAL DIAGNOSTIC

WWW.GL-BIOCONTROL.COM