

PROTOCOL

REHYDRATION OF PHAGE Φ X174

Description


Reference: PHIX174

Biosafety Level: 1

Storage temperature: < -20°C

For research use only.

Rehydration of the phage Φ X174



You get a ready-to-use viral suspension.



After use, close the aluminium bag to avoid rehydration of the freeze-dried viruses. Store in the freezer (< -20°C).

Growth conditions

Method: ISO 10705-2

Analysed volume: 5 ml

Media: MSA

Incubation temperature: 36 ± 2°C

Incubation time: 18 ± 2h

Atmosphere: aerobic