



DENDRIDIAG® IW

ATP tests for industrial water



MEASURE THE QUANTITY OF BACTERIA IN YOUR INDUSTRIAL WATER CIRCUIT IN UNDER 2 MIN

SAVE

up to 40% on treatment products. Avoid financial penalties and production shutdowns.

SECURE

your water circuit by limiting the risks associated with *Legionella* and by anticipating microbial shifts.

COMPLY

with the regulations. Have a reliable microbial indicator at your disposal.

MASTER YOUR WATER CIRCUIT



Monitor

the microbial quality of your water



Identify

the critical points in no time at all



Take action

immediately in a crisis situation



Optimize

your cleaning and disinfection process

THEY TRUST US



Pierre B. - **ENGIE**

I have been using ATP tests for the monitoring of cooling towers and domestic hot water systems for 4 years now. **ATP tests have become essential for us.** It is a great tool for use on site; very easy-to-use and rapid to anticipate shifts and identify critical areas.



Laurence F. - **BWT**

I use the DENDRIDIAG® kit to demonstrate the effectiveness of my biocide products. In some cases, I also use it **to improve my treatment strategies**, which reassure my clients. It also allows me to meet cooling tower regulations.

For more information, please visit ► www.gl-biocontrol.com/en/IW




WHAT IS DENDRIDIAG®?



A complete kit to measure the quantity of bacteria in industrial water


◀ Watch the video

- Cooling water systems,
- Industrial processes,
- Electroplating baths...




Get the result in a second thanks to the web app!

A FAST AND RELIABLE ANALYSIS



INDIVIDUAL SUPPORT




KIT CONTENT



The reagents

- 60 measurement tubes LUMITUBE
- 1 EXTRACTANT dropper
- 1 STANDARD dropper



The consumables

- 60 sampling containers
- 60 10 ml syringes
- 60 0.45 µm filters



The measuring device

- 1 luminometer KIKKOMAN SMART
- 1 microtube rack
- 1 USB flash drive

NEED INFORMATION? CONTACT US



Head office

5, avenue de l'Europe - Hélioparc
34 830 CLAPIERS (FRANCE)

Grand Est Agency

8, rue Alfred Kastler
Maison de l'entreprise
54 320 MAXÉVILLE (FRANCE)



+ 33 (0)9 67 39 35 20



contact@gl-biocontrol.com



Join us on LinkedIn

GLBIOCONTROL

RAPID MICROBIAL DIAGNOSTIC

WWW.GL-BIOCONTROL.COM